

# ELVeS® Radial® 2ring swift



## Why to use ELVeS® Radial® Swift Fibers?

- Just 5Fr Introducer Size
- or Easy access via 14GA catheter
- Well known FUSION® safety concept

The swift fibers are designed to complete our product range. You decide upon access by direct puncture via 14GA canula or using a ordinary 5FR Introducer set (instead of 6Fr).


Adapt flexible to the particular vein in terms of energy, power and approach.



Placing the catheter



Entering the vein



Starting laser treatment

### Fibers\*

REF	Product	PU **	length [m]	ø fiber tip [mm]	CI***
503100140	ELVeS® Radial® swift, IC	10	2.5	1.5	5FR/14GA
503100142	ELVeS® Radial® swift 2ring, IC	10	2.5	1.5	5FR/14GA

### Kits\*

503100141	ELVeS® Radial® swift Kit, IC	10	2.5	1.5	5FR/14GA
503100143	ELVeS® Radial® swift 2ring Kit, IC	10	2.5	1.5	5FR/14GA

\* Only to be used with 1470 nm \*\* Packaging unit \*\*\* Compatible introducer

**Imprint**  
biolitec AG  
Untere Viaduktgasse 6/9  
A-1030 Wien  
Phone: +43 1 3619 909 50  
info@biolitec.de  
www.biolitec.com

All fibers are free of latex and DEHP.  
Our fibers are single use products  
(unless otherwise indicated) delivered  
sterile for immediate use.  
biolitec®, LEONARDO®, FUSION® and  
ELVeS® are registered trademarks  
owned by biolitec.

Contact us



**bio**  
**LITEC**®